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
The screenshot shows the Instadose user interface. At the top, it says 'Inicio > Panel de información' and 'You are logged in as Adrian Milei Log Out'. The main header says 'Bienvenido Adrian'. Below this is a 'Visualización rápida' section with a table of user information:

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

The screenshot shows the main dashboard of the instadose application. The user is logged in as Adrian Milei. The dashboard includes a navigation menu on the left with options for HOME, SUPPORT, and MY ACCOUNT. The main content area is titled 'Bienvenido Adrian' and features a 'Visualización rápida' section with a table of user and device information. A 'Leer dispositivo' button is highlighted with a red box.

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
The screenshot shows the instadose web application during a device reading process. The user is logged in as Adrian Milei. The main content area displays a progress bar and a message: 'Espere hasta que se envíen los datos al servidor...'. A 'Leer dispositivo' button is visible. A window titled 'Instadose(tm) Device Reader Ver 1 Rel 2.04' is open, showing a progress bar and a message: 'Instadose Reader Devices Detected=1, Maximum=999'. The window also contains several icons representing different functions.

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You are logged in as Adrian Milei [Log Out](#)

HOME  
SUPPORT  
**MY ACCOUNT**

Service provided by:  
**MIRION** Dosimetry Services  
TECHNOLOGIES Division

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### Radiation in Daily-life

Radiation Dose  
(Millisieverts: mSv per year)

250.00	Upper limit of radiation dose permitted for people who engage in emergency work
50.00	Upper limit of radiation dose permitted for radiation workers, police, and firefighters who engage in disaster prevention
10.00	Radiation dose in Guarapari (Brazil) per year
10.00	Abdomen or pelvic CT scan (10.00 mSv/ each time)
6.90	Chest CT scan (6.90 mSv/ each time)
2.40	Natural occurring radiation global average received from space, earth, air and ingestion
1.00	Dose limit for public per year (except for medical care)
0.60	Gastrointestinal x-ray examination (60 mrem/ each time)
0.20	Round-trip flight between Tokyo and New York
0.13	Mammogram x-ray examination (0.13 mSv/ each time)
0.06	Limb and joint x-ray examination (0.06 mSv/ each time)
0.05	Chest x-ray examination (0.05 mSv/ each time)
0.05	Standard dose of radiation around a nuclear plant light waste reactor)

[Legal](#) | [Terms and Conditions](#)

Copyright © 2014 Mirion Technologies, Inc. All Rights Reserved.

## Nuclearlab SRL

Dr. Luis Beláustegui 776 (1416) CABA, Argentina Tel – Fax: +54-11-4584-8308  
[instadose@nuclearlab.com](mailto:instadose@nuclearlab.com) [www.nuclearlab.com](http://www.nuclearlab.com)